

EARTH 48+48 TWIN

Aluminum Body

LED

Number of LED
48 + 48



Intensity
160,000 x 2 Lux



Color Temp
3800 K - 5200 K



Power
48 + 48 Watt



Diameter Of Light
500 mm + 500 mm



Multi functional Controller
With Remote

Technical Data

Number of LED	48 + 48	Focusing	Adjustable
Intensity	160, 000 x 2 Lux	LED Average Life	> 50000 Hrs.
Spot Diameter	100 - 150 mm	Intensity Control	10 -100 Digital
Color Temperature	3800 - 5200K	Power Supply	220V/50Hz AC
Depth illumination	10-15 cm	Power	48 + 48 Watt
Color Rendering Index	95 Ra	Optional	Battery Backup
Diameter of Light	500 + 500 mm		360 degree structure
Special Function	5 Memory Function with Remote		

Special Features

- Perfect Compatibility of Two Light heads.
- Easy Positioning with Detachable Sterilized Central Control Handle
- Fast and Simple Operation with Remote or Control Panel